## Caring for you

## Contact Us:





Neusoft Medical Systems



in Neusoft Medical



## **NeuAngio-CT**

One-stop Multi-modality Solution





NeuAngio-CT, a pioneering one-stop multi-modality solution for interventional diagnosis and treatment, takes "fast, precise, safe and flexible" as the design concept, effectively shortens the treatment path and seamlessly links diagnosis and treatment through the synergy of high-end sliding CT and DSA.

NeuAngio-CT helps the users to save valuable time, optimize procedure efficiency and accuracy. It also avoids the risk of patient transfer, reduces the radiation dose received by the operators and the patients, and improves patient outcomes.



\*1 Table-side Touch Screen shown in the brochure is optional \*2 58" large screen shown in the brochure is optional

Neusoft Medical Systems Co., Etd.





## ▼ Railless 7-axis Ceiling Mounted Angio System




- · 65 °/s rotational angiography saves more time for the treatment
- · Unique lateral movement is beneficial to radial artery interventional procedures

## ▼ Fast Sliding CT with High Precision



- · Advanced 128-slice Sliding CT
- · 200mm/s sliding speed and 0.374s scan speed to achieve faster image acquisition



- · Fast switching and fusion to obtain more accurate image information
- · Both systems enables remote parking and no interfere with each other
- · Integrated user interface, displaying full image of the operation

Neusoft Medical Systems Co., Ltd.

One-stop Multi-modality Solution



### Glorious Atlas ▼




Transparent Vessel ▼

Dual-Volume Rendering ▼


Tumor Evaluation ▼




Glorious Dual Vessel ▼




Glorious CT ▼





Body Perfusion ▼



Lung Nodules ▼




3D Reconstruction ▼ Head and Neck CTA ▼



Multi-Volume Rendering▼ MPR▼



Lung Density ▼



Neusoft Medical Systems Co., Ltd. 05/06

Tumor Intervention Solution

Cardiac Intervention Solution

Head Plain Scan ▼

Brain Perfusion ▼



Lower Limb CTA ▼




Coronary Tree VR\*▼

Aortic CTA ▼







TAVR\* ▼






Cardiac Calcium Analysis\* ▼



Cardiac Function Analysis\* ▼



Coronary Analysis\* ▼

TRO\*▼

Right Coronary CPR\* ▼

Aortic Dissection ▼



One-click Cardiac Segmentation\*▼



\* Optional only

Neusoft Medical Systems Co., Ltd.



**D-MultiPulse** 

\* Optional only


中

D-5 Elements AEC

# Flexible

Sliding CT is designed to be used both in CT Examination Room and NeuAngio-CT operating room to maximize the utilization.











**NeuAngio-CT Main Console** 



NeuAngio-CT One-stop
Post-processing Workstation

# Flexible



Fast switching and fusion of Sliding CT/DSA to save more time for patients in emergency



Maximizing of CT utilization rate and optimizing return of investment of hospitals



Customized installation is provided to release larger table-side operating space and more flexible positioning



## Neusoft Global Service & Logistics Network



### After-sales service and support

- Remote service capabilities bring Neusoft expertise to you IMMEDIATELY, no matter where you are!
- Identifying and correcting PROMPTLY and PROACTIVELY, minimizing downtime and patient inconvenience.
- Global logistics network enables fast response regarding parts and supplies.
- \* Note: The contents of this publication and the listed parameters are for reference only and not intended as legal offers or commitments. Neusoft Medical Systems reserves the right to modify the contents, design, specifications and options described herein without prior notice, and will not be liable for any consequences resulting from the use of this publication. Please contact your local Neusoft sales representative for the current information. The specific sales product configuration is subject to the actual contract signed by Neusoft.

\* Not available in the United States.

Neusoft Medical Systems Co., Ltd.